

Veratox®

For DON 5/5



Intended Use

Veratox for DON 5/5 is used for the quantitative analysis of DON in grain and grain products, such as wheat, wheat flour, wheat midds, wheat bran, corn, cornmeal, corn screenings, barley, malted barley and oats.

The Test

The test is a competitive direct ELISA that provides exact concentrations in parts per million (ppm). Free DON in the samples and controls competes with enzyme-labeled DON (conjugate) for the antibody binding sites. After a wash step, substrate reacts with the bound conjugate to produce blue color. A microwell reader is used to yield optical densities. Control optical densities form the standard curve, and sample optical densities are plotted against the curve to calculate the exact concentration of DON.


The Procedure

Samples must be extracted prior to testing.

1. Add 100 μ L conjugate to each mixing well.
2. Add 100 μ L controls and samples to mixing wells.
3. Mix. Transfer 100 μ L to antibody wells. Incubate for 5 minutes.
4. Dump liquid from antibody wells.
5. Wash wells thoroughly with deionized water.
6. Tap out water on paper towel.
7. Transfer 100 μ L substrate from reagent boat to antibody wells using 12-channel pipettor. Incubate for 5 minutes.
8. Transfer 100 μ L Red Stop from reagent boat to antibody wells.
9. Read results using a microwell reader with a 650 nm filter.

See package insert for complete instructions

Results



Product Specifications

Lower limit of detection:	0.1 ppm
Range of quantitation:	0.25 – 2.5 ppm
Controls provided:	0, 0.25, 0.5, 1, 3 ppm
Testing time:	10 minutes
Antibody cross-reactivity:	DON 100%. Also cross-reactive with 3-acetyl DON but does not bind to other trichothecenes.
Tests per kit:	Up to 38
Approvals:	USDA/GIPSA #2002-106

Materials Recommended But Not Provided

Available from Neogen*

1. Graduated cylinder
2. Sample collection cup, or container with 250 mL capacity
3. Neogen filter syringe, Whatman #1 filter paper, or equivalent
4. Sample collection tubes
5. Agri-Grind grinder, or equivalent
6. Scale capable of weighing 10-50 grams
7. Microwell reader with a 650 nm filter
8. 100 μ L pipettor
9. 12-channel pipettor
10. Tips for 100 μ L and 12-channel pipettors
11. Microwell holder
12. Timer
13. Wash bottle
14. 2 reagent boats for 12-channel pipettor

Not Available from Neogen

1. Paper towels or equivalent absorbent material
2. Plastic $\frac{1}{2}$ gallon bucket for use as waste receptacle
3. Distilled or deionized water
4. Waterproof marker

* See *Equipment and Accessories*, pages 72-74

Ordering Information

Prod#.	Product description
8331	Veratox for DON 5/5
8052	Mycotoxin Extraction Kit