


S/J PICKING TRAY

The S/J Picking Tray was designed to make the job of a grain grader easier by turning all kernels over at one time. The grain is placed on the corrugated half of the tray and with a few short, quick shakes the grain is spread in one layer. With the use of tweezers, kernels can be "plucked" from the formed rows without disturbing the kernels on either side of the one being removed. After all rows have been reviewed, the second half of the picking tray (foam side) is now placed on top of the first half (corrugated half). Holding the two trays tightly together, the user turns the trays over and now places the foam side of the tray on the bottom. The kernels remain lined-up and the opposite side of the grain can now be viewed for further grading. SIB-L Ship wt. 6 lbs., Dims 24"L x 19"W x 4"H. SIB-S Ship wt. 3 lbs., Dims 17"L x 12"W x 3"H. SGT Ship wt. 4 lbs., Dims 24"L x 7"W x 7"H.

MADE IN USA

Seedburo Inspection Board - Large (corn, soybean, sunflower)	No. SIB-L
Seedburo Inspection Board - Small (barley, wheat, lentils)	No. SIB-S
Seedburo grain trough	No. SGT



SEEDBURO FORCEPS

Our forceps are used for analyzing grain or seed and can be used in conjunction with our inspection boards as featured above. Our medium serrated points are suggested for use with larger seed and the fine points for more delicate seed.

Nickel Plated, 5-1/4" medium sharp point	No. 661
Nickel Plated, 6" blunt point	No. 662
Stainless Steel, 4-3/4" sharp point	No. 663
Stainless Steel, 5-1/2" sharp point with curved tip	No. 664
Nickel Plated, 8" blunt point	No. 665
Nickel Plated, 6-1/2" sharp serrated point, cross locking	No. 667
Nickel Plated, 6" sharp serrated point with curved tip	No. 668
Chromium Plated, 6-3/8" serrated medium sharp point	No. 669
Nickel Plated, 10" blunt point with curved tip	No. 670
Nickel Plated, 10" blunt point	No. 671

FIVE PIECE ALL-PURPOSE TWEEZERS KIT

A complete, basic tweezers set. An ideal way to sample and select the proper type tweezer for most applications. Multiple tips available for the task at hand. Includes deluxe vinyl pouch for handy storage.



Five Piece All-Purpose Tweezers Kit	No. 672
---	----------------

SEEDBURO RAZOR BLADES

These stainless steel, single edge, 3 facet razor blades are a "must have" for any seed technologist. Features an aluminum back, PTFE (Teflon) coating and are individually protected. Length: 1.56" (39.6mm), Cutting Edge Length: 1.527" (38.8mm), Width: 0.77" (19.6mm), Thickness: .009" (.220mm). Sold in box of 100 blades. Wt. 1 lb.

Single Edge, PTFE Coated Stainless Steel Blades, 100/box....	No. 121-3
--	------------------

MEDIUM & HEAVY DUTY BENCH KNIVES

A handy tool for cutting seeds, digging under seed coats and many other uses. These bench knives are available in medium and heavy duty (two sizes). Blade is firmly set in hardwood handle.



#39110

Medium Duty Bench Knife, 1-3/4" slim blade, 5" overall	No. 39180
Heavy Duty Bench Knife, 1-3/4" blade, 5-3/4" overall	No. 39105
Heavy Duty Bench Knife, 2-5/8" blade, 6-5/8" overall	No. 39110