


SEEDBURO SAMPLE CONTAINERS AND PANS



QUART AND PINT SAMPLE BOTTLES

Seedburo's Quart and Pint Sample Bottles are made of a durable polyethylene. Includes an enameled plastic screw-top cap that has a gasket insert to provide an air-tight seal. The No. 630 Quart Bottle meets the specifications of the USDA/ GIPSA (FGIS) for official inspection to minimize moisture loss and preserve quality while in short term storage. Both units feature a 1-3/4" neck opening.

- Quart Sample Bottle, 3-3/4" L x 3-3/4" W x 7" H **No. 630**
- Pint Sample Bottle, 2-3/4" Dia x 6" H **No. 635**



HALF GALLON AND GALLON CONTAINERS

Our Half Gallon and Gallon Containers are made of a durable, high density polyethylene. The 5095 Half Gallon Container features an air-tight, snap-on lid. Commonly used in inspection offices and laboratories, they are great for "stackable" storage. Both the 625 and 640 containers have screw-top plastic cap. **MADE IN USA**

- Half Gallon Square Container, Snap-On Lid, 5" L x 5" W x 6-3/4" H **No. 5095**
- Half Gallon Container, Screw-Top Cap, 5" L x 5" W x 6-1/2" H **No. 625**
- Half Gallon Round Container with Snap-on Lid, 6" dia. X 5-3/4" H **No. 645**
- Gallon Container, 6" Dia x 10" H with 3-3/4" neck opening, Screw-Top Cap **No. 640**



TRIANGULAR SAMPLE PANS

The No. 64 Sample Pan is made of mild steel and has an enameled, dark blue color. This pan has become a standard and features a tapered pour spout. The No. 64P is the same pan made of durable, white, static resistant plastic. **MADE IN USA**

- Blue Tin Triangular Sample Pan, 10" L x 10" W x 2-1/2" H **No. 64**
- White Plastic Triangular Sample Pan, 10" L x 10" W x 2-1/2" H **No. 64P**



SPOUTED SAMPLE PANS

These pans are made of aluminum and are ideal for controlled pouring of samples. The No. 33 and No. 31 feature a funnel-type pour spout. **MADE IN USA**


- Spouted Sample Pan, 1-1/2 Quart 11-3/8"L x 7-1/4" W x 1-3/8"H **No. 33**
- Spouted Sample Pan, 3 Quart 15-1/4"L x 9-7/8"W x 2-1/8"H **No. 31**
- Mini-Tin Pan (7" L x 5" W x 7/8" H) **No. 7**



ROUND SAMPLE PANS

The 44 sample pan series is made of stainless steel and are available in the following sizes:

- Round Aluminum Sample Pan (Pkg/12) - 5" Dia x 1" H..... **No. 44**
- Round Sample Pan, 2 qt liquid (1.7 qt dry) - 8" Dia x 2-3/4" H **No. 44-2**
- Round Sample Pan, 4 qt liquid (3.4 qt dry) - 9-1/4" Dia x 2-7/8" H **No. 44-4**
- Round Sample Pan, 6 qt liquid (5.1 qt dry) - 11" Dia x 3-1/2" H **No. 44-6**



SEEDBURO ANTI-STATIC WEIGH BOATS

These durable, lightweight weigh boats are ideal for use with digital balances or small weighing applications. Resist temperatures to 200°F and flex for easy pouring. Rubberized polystyrene with anti-static properties. Designed with shallow shape and flat bottom to resist tipping. Sold in packages of 50 each. Also available in blue.

- Small Weigh Boat (1-3/4" L x 1-3/4" W x 3/8" D), White **No. 1420S**
- Medium Weigh Boat (3-1/4" L x 3-1/4" W x 1" D), White **No. 1420M**
- Large Weigh Boat (5-1/2" L x 5-1/2" W x 7/8" D), White **No. 1420L**