



CIRCLINE MAGNIFIER LAMP

Seedburo's Circline Magnifier Lamp provides uniform lighting to illuminate and enlarge objects placed under the lens. Unit is ideal for laboratories, quality control and grain inspection. Built-in handle lets you maneuver the lamp easily. The 5-inch, 3-diopter lens increases magnification 75% at a focal length of 13". The arm reaches a full 28". Circular lamp uses a 22-watt bulb. Clamp-On Model Ship wt. 12 lbs, Ship dims. 25" L x 24" W x 5" H. Desk Base Model Ship wt. 28 lbs, Ship dims. 13" L x 13" W x 25" H. **MADE IN USA**

Circline Magnifier:

- Clamp-On Model 110V, 60 Hz.....**No. MC150**
- Clamp-On Model 220V, 50Hz**No. MC150/C**
- Desk Base Model 110V, 60Hz.....**No. MC100**
- Desk Base Model 220V, 50Hz**No. MC100/C**
- Replacement Bulb, 22 Watt (for 110V and 220V) **No. MCB**
- Add-A-Lens attachment to increase magnification to +275% at a focal length of 3-3/4" **No. M1400**



HIGHLIGHTING MAGNIFIER LAMP

The Highlighting Magnifier is perfect for highlighting flaws, imperfections and other irregularities on the object being inspected. Magnifier features a 5-inch, 3-diopter lens to increase magnification to 75% at a focal length of 13". Arm reach is 26" and the lamp uses a 13-watt fluorescent bulb (9-watt for 220V). Clamp-On Model Ship wt. 12 lbs, Ship dims 25" L x 24" W x 5" H. Desk Base Model Ship wt. 28 lbs, Ship dims. 13" L x 13" W x 25" H. **MADE IN USA**

Highlighting Magnifier:

- Clamp-On Model 110V, 60Hz..... **No. MG450**
- Clamp-On Model 220V, 50Hz**No. MG450/C**
- Desk Base Model 110V, 60Hz..... **No. MG400**
- Desk Base Model 220V, 50Hz**No. MG400/C**
- Replacement Bulb:**
- 13 Watt (for 110V) **No. MGB**
- 9 Watt (for 220V) **No. MGB/C**



SWIFT STEREO MICROSCOPES

These modular stereo microscopes offered by Seedburo incorporate the latest optical technology. The SBM80 Series is the best choice for most microscopic applications. Features include: 360° rotating body, two selectable powers, super-wide 10X eyepieces, wide interpupillary range, diopter equipped eyetube and Ultra-Lite illuminating system. 110V, 50/60Hz. Ship wt. 15 lbs., Ship dims. 21" L x 21" W x 16" H.

Model	Objective	Eyepiece	Total Magnification
SBM91	1X, 2X	10X	10X & 20X
SBM95	1X, 3X	10X	10X & 30X
SBM98	2X, 4X	10X	20X & 40X

SEEDBURO HAND-HELD MICROSCOPE

The IPM Scope combines a digital movie camera, precision optics and LED lighting in a microscope that fits in the palm of your hand. 40 to 140X magnification lets you zoom in on the fine details of insects, textiles, electronics or anything else too small to see.

Place the IPM Scope over the object and comfortably view the image on your computer screen, instead of straining to look into a tiny eyepiece. Capture images of movies, easily add labels, make measurements and even draw right on the live video image. The IPM Scope connects to your computer's USB port where it's always ready to use. Ship wt. 2 lbs, Ship dims. 7"x5"x3"

- Hand Held Microscope & Digital Camera with Imaging Software..... **No. IPM-SCOPE**
- Scope Stand..... **No. IPM-STAND**