

SEEDBURO PC400 SERIES TOP-LOADING BALANCE

The Seedburo PC400 Top-loading Balances have the exclusive quad-spring design that virtually eliminates damage from shock overload. The scales read in pounds, ounces, kilograms or grams. They run on an internal rechargeable battery or AC plug-in power adapter. The scale can even be used while recharging. The small size of the PC400 allows you to use it in tight spaces. In addition, the scales are easy to clean and made of 304 stainless steel, which is wipe down safe. Complete with bi-directional RS232 output. The PC400-05 and PC400-10 are both NTEP approved*. Power requirements 115V or 230V, 50/60Hz. Actual dims 8.75" L x 8.5" W x 4.25" H.



PC400 Balance

Model	Capacity x Readability		
	Grams	Ounces	Pounds
No. PC400-02	1000 x 0.1	35 x 0.005	2 x 0.0002
No. PC400-05*	2300 x 1	80 x 0.05	5 x 0.002
No. PC400-10*	4500 x 2	160 x 0.1	10 x 0.005
No. PC400-Adapter	Adapter for 220V Operation of above scales.		

OHAUS COMPACT SCALES

Ohaus Compact Scales are ideal when cost matters and dependable results are a must. These uniquely designed scales have a wide variety of applications, more than any other scale in their class. The three portable models range in capacities from 200g to 5000g, and are perfect for use in light production, in the field, the kitchen, the office and more. Easy operation, long battery life, and great performance make Ohaus Compact Scales a truly remarkable value. Standard features are a low battery indicator, overload and under load indicators, auto shut-off, rugged plastic construction and tactile feedback keypad. Weighing modes include grams, pounds, and ounces. The CS series is powered by 3-AA batteries (included) or an optional AC adapter (120V or 220V, 50/60Hz), Ship wt. 2 lbs. Ship Dims 9½" L x 7" W x 2¼" H.



CS200 Scale

Model	Capacity	Readability	Platform
No. CS200	200g	0.1g	7.1" Diam.
No. CS2000	2000g	1.0g	7.1" Diam.
No. CS5000	5000g	1.0g	7.1" Diam.
No. 90524-66	110V, 50/60Hz Adaptor		
No. 90524-63	230V, 50/60Hz Adaptor (European Type)		

OHAUS ADVENTURER™ PRO PRECISION ELECTRONIC BALANCES

The new Ohaus Adventurer™ Pro continues to combine unsurpassed value with solid reliability. Precision performance, versatility, and cutting edge desing come together to deliver accurate results. Adventurer™ Pro's unique design intergrates a custom two-line backlit LCD for easy viewing and simple menu navigation. An ergonomic leveling system with an upfront bubble indicator quickly ensures that the balance is level. **Standard Features/Equipment:** Large stainless steel weighing platform, two-line backlit LCD, bi-directional RS232 interface, leveling system, in-use cover, intergal security bracket, weigh below hook, GLP data output options, adjustable environmental settings, stability indicator, overload and underload indicators, low battery indicator, auto shut-off, AC adapter, draftshield on Model AV53.



Adventurer Pro

Model No.	AV212	AV412	AV812	AV811	AV2101	AV4101	AV8101	AV8100	
Capacity (g)	210	410	810	810	2100	4100	8100	8100	
Readability	0.01			0.1				1	
Weighing Units	Milligram, Gram, Kilogram, Ounce, Pound, Carat, Pennyweight, Ounce Troy, Grain, Newton, Hong Kong Tael, Singapore Tael, Taiwan Tael, Momme, Tical, Baht, Mesghal, Tola, Custom								
Tare Range	To Capacity by Subtraction								
Stabilization Time (s)	1.5	3			2				
Operating Temp Range	50° - 86° F/10° - 30° C								
Power Requirements	AC Adapter (Included) or 4 AA Batteries (Not Included)								
Calibration	Digital with External Weight								
Display Type	2-Line LCD w/Backlight								
Pan Size (in/cm)	4.7 / 12 dia.		4 x 1 / 10 x 2.5						
Dimensions (in/cm)	7.5 W x 2.8 H x 10 D / 19.3 W x 7.2 H x 25.4 D								
Net Weight (lb./Kg)	4.8 / 2.2		5.2 / 2.4			5.5 / 2.5			