

## COUNT-A-PAK® SEED TOTALIZER

Counts Seeds in a Given Weight or Volume Quickly, Accurately and Simply.

This counter has the potential to save producers valuable seed, time, labor and money in ensuring bagging scales are properly calibrated each day. Determine seed count per pound or 1000 kernel weights with the 77 Totalizer. The speed of the feeder is adjusted by a built-in control on the feeder base. With 2" LCD readout, the counter displays total count of the seed which is easily read under most lighting conditions. Available with a 7", 10" or 10" shallow bowl, you'll find a suitable counter for almost any need.



Power: 175 Watts; Actual Dims Control Box: 15" L x 10" W x 5" H; Actual Dims Feeder Base: 9" L x 9" W x 9" H; Net Wt. (both pieces) 32 lbs.; Ship Dims (complete): 32" L x 16" W x 14" H; Ship Wt. (complete): 55lbs.

Seed Totalizer with 10" bowl 110V/60Hz .....	No. 77-10/B
Seed Totalizer with 10" bowl 110V/50Hz .....	No. 77-10/D
Seed Totalizer with 10" shallow bowl 110V/60Hz .....	No. 77-10S/B
Seed Totalizer with 10" shallow bowl 110V/50Hz .....	No. 77-10S/D
Seed Totalizer with 7" bowl 110V/60Hz .....	No. 77-7/B
Seed Totalizer with 7" bowl 110V/50Hz .....	No. 77-7/D
Transformer for 220 volt operation .....	No. K25
Spare parts kit for 10" bowl 110 or 220V .....	No. 77-10SPK*
Spare parts kit for 7" bowl 110 or 220V .....	No. 77-7SPK*
Replacement Frosted Glass Only (each) .....	No. 77FG
Replacement Clear Glass Only (each) .....	No. 77CG
Domestic Packing .....	No. 77/DOM

\* Kit consists of replacement springs, frosted and clear glass.

### Quick Disconnect Feeder Bowl Standard on Model 801 only.



Unlock Mechanism...



Remove and Dump...



Replace and Lock

The following chart gives approx. bowl capacities, average feed rates\* during batch counting process, and time required to empty a full bowl.

Variety	Average Speed	10-Inch Capacity	Time to Empty	7-Inch Capacity	Time to Empty	10S-Inch Capacity	Time to Empty
Corn	160	3,600	20-25 min	650	4-5 min	1,500	9-11 min
Soybeans	200	7,000	30-40 min	1,200	5-7 min	2,300	12-15 min
Popcorn	225	9,000	35-45 min	1,500	6-8 min	3,500	13-17 min
Wheat	265	30,000	1 ¾-2 hrs	4,500	15-20 min	10,000	35-45 min
Milo	315	40,000	2-2 ¼ hrs	5,000	15-20 min	11,500	39-45 min

\* The average speed (seed/minute) shown are for batch counting. Totalize counting can be 2 to 3 times faster. Capacity is approximate number of seeds in a full bowl.

**U.S. / Canada**  
Toll Free:

**1-800-284-5779**

**Fax: 312-738-5329**

**International: 312-738-3700**