



CIRCLINE MAGNIFIER LAMP

Seedburo's Circline Magnifier Lamp provides uniform lighting to illuminate and enlarge objects placed under the lens. Unit is ideal for laboratories, quality control and grain inspection. Built-in handle lets you maneuver the lamp easily. The 5-inch, 3-diopter lens increases magnification 75% at a focal length of 13". The arm reaches a full 28". Circular lamp uses a 22-watt bulb. Clamp-On Model ship wt 16 lbs, dims 13" L x 12" W x 34" H. Desk Base Model ship wt 28 lbs, dims 25" L x 13" W x 17" H.

Circline Magnifier:

110V, 60 Hz, Clamp-On Model	No. MC150
220V, 50Hz, Clamp-On Model	No. MC150/C
110V, 60Hz, Desk Base Model	No. MC100
220V, 50Hz, Desk Base Model	No. MC100/C
Replacement Bulb, 22 Watt (for 110V and 220V)	No. MCB
Add-A-Lens attachment to increase magnification to +275% at a focal length of 3¾"	No. M1400



HIGHLIGHTING MAGNIFIER LAMP

The Highlighting Magnifier is perfect for highlighting flaws, imperfections and other irregularities on the object being inspected. Magnifier features a 5-inch, 3-diopter lens to increase magnification to 75% at a focal length of 13". Arm reach is 26" and the lamp uses a 13-watt fluorescent bulb (9-watt for 220V). Clamp-On Model ship wt 15 lbs, dims 13" L x 12" W x 34" H. Desk Base Model ship wt 27 lbs, dims 13" L x 13" W x 25" H.

Highlighting Magnifier:

110V, 60Hz, Clamp-On Model	No. MG450
220V, 50Hz, Clamp-On Model	No. MG450/C
110V, 60Hz, Desk Base Model	No. MG400
220V, 50Hz, Desk Base Model	No. MG400/C

Replacement Bulb:

13 Watt (for 110V)	No. MGB
9 Watt (for 220V)	No. MGB/C

SEEDBURO PORTABLE DUAL-PURPOSE SCOPES

These Dual-Purpose Scopes are lightweight and portable. Constructed of attractive and durable aluminum and plastic. Well suited for laboratory and field use, these scopes are ideal for macro- and microscopic quality control and inspection. Use as a microscope for grain inspection, identifying insects or to verify coating procedures on seeds. Use as a monocular for distance viewing. Each scope includes 3" tall micro lens stand, neck strap, lens cover, and hard carrying case with strap. DPS20 dims (without stand) 2.9" H x 1.3" dia, ship wt .5 lbs. DPS40 dims (without stand) 5.1" H x 1.7" dia, ship wt 2 lbs.



As a Microscope:

Model No.	Magnification (Highest)	Field of View (Max)	Working Distance
DPS 20	20 X	13.0mm	3.1"
DPS 40	40 X	8.0mm	2.5"

As a Telescope:

Model No.	Magnification	Field of View	Close Focus Distance	Objective Lens Aperture
DPS 20	6 X	10°	10.2"	16mm
DPS 40	8 X	8.5°	23.6"	30mm

SWIFT STEREO MICROSCOPES

These modular stereo microscopes offered by Seedburo incorporate the latest optical technology. The SBM80 Series is the best choice for most microscopic applications. Features include: 360° rotating body, two selectable powers, super-wide 10X eyepieces, wide interpupillary range, diopter equipped eyetube and Ultra-Lite illuminating system. Options available for greater magnification. 110V, 50/60Hz. Ship wt 13 lbs., dims 21" L x 15" W x 11" H.



Model No.	Objective	Eyepiece	Total Magnification
SBM91	1X, 2X	10X	10X & 20X
SBM95	1X, 3X	10X	10X & 30X
SBM98	2X, 4X	10X	20X & 40X

U.S. / Canada
Toll Free:

1-800-284-5779

Fax: 312-738-5329

International: 312-738-3700