



Bunn® Grinding Mill

The Bunn Grinder is used for the preparation of corn and other grain samples for procedures such as Mycotoxin testing. The large capacity hopper holds up to 3 lbs of Corn for efficient grinding of large sample lots. The turbo action draws the sample through the grinding chamber for greater grinding efficiency and will process approximately two pounds of corn per minute. Selectable grind range from coarse to fine. A convenient cleaning lever ensures complete clean out from the grinding chamber. The Bunn Grinder is built with a dependable ¾HP motor which powers the heavy duty precision fine burrs for a quality grind that cuts the kernel rather than crushes it. Burrs carry a 4-year warranty. Unit draws 9 amps and therefore power requirements are for a minimum 15 amp circuit. Net wt: 50 lbs, Actual dims: 17" (L) x 8" (W) x 28" (H), Ship wt: 56 lbs, Ship dims: 23" (L) x 14" (W) x 33" (H).

- Bunn Grinding Mill, HD Black Commercial, 120V, 60Hz..... **No. 22100**
- Burr Set Kit, Complete **No. 5861**

22100 with Corn sample

Seedbuero Electric Grinding Mill

Seedbuero's Electric Grinding Mill offers a fast means for grinding your product. The unit is powered by a 1/3HP motor with an output speed of 90 RPM. Ideal for the laboratory or mill house, this grinder features the same type of "grinding head" as the Seedbuero hand grinder mounted on a high speed motor for faster and more consistent results. Output capacity of most products is 40 lbs/hour. The maximum fineness depending on product and moisture content: 100 mesh. Hopper size: 4½" W x 5" D. Counter top to discharge height: 6". (230V model with gear box drive).



115V Unit: Net wt: 55 lbs, Actual dims: 25" (L) x 8" (W) x 12" (H). Ship wt/dims: **Carton 1 (Base):** 51 lbs, 20" (L) x 12" (W) x 12" (H). **Carton 2 (Head):** 12 lbs, 14" (L) x 11" (W) x 11" (H).



230V Unit: Net wt: 65 lbs, Actual dims: 28" (L) x 8" (W) x 12" (H). Ship wt/dims: **Carton 1 (Base):** 71 lbs, 24" (L) x 12" (W) x 12" (H). **Carton 2 (Head):** 12 lbs, 14" (L) x 11" (W) x 11" (H). Also Available in Stainless Steel.

- Grinding Mill with Both Feeds, 115V, 60Hz.....**No. 438B/B**
- Grinding Mill with Tooth Feeds, 115V, 60Hz.....**No. 438T/B**
- Grinding Mill with Worm Feeds, 115V, 60Hz..... **No. 438W/B**
- Grinding Mill with Both Feeds, 230V, 50Hz.....**No. 438B/C**
- Grinding Mill with Tooth Feeds, 230V, 50Hz.....**No. 438T/C**
- Grinding Mill with Worm Feeds, 230V, 50Hz..... **No. 438W/C**



Seedbuero Hand Grinder

The Seedbuero Hand Grinder features an adjusting screw to set the grinder to a wide range of fineness. The unit is made of high grade cast iron with an approximate capacity hopper of 2 lbs. Interchangeable feeds available. Tooth feed for dry grinding or worm feed for wet and/or oily grinding. This grinder is ideal for preparing soft or hard grains for cereal and bread, grinding Corn for Aflatoxin testing or for general purpose grinding on a wide variety of commodities such as commercial grains, Coffee beans, Peanuts and other nuts, meats, etc. Seedbuero Hand Grinders are used by government experiment stations, research laboratories and many grain elevators, feed mills and flour mills. Net wt: 12.5 lbs, Ship wt: 16 lbs, Actual dims: 14" (L) x 7" (W) x 15" (H), Ship dims: 13" (L) x 13" (W) x 13" (H).

- Hand Grinder with Both Feeds.....**No. 438HB**
- Hand Grinder with Tooth Feeds.....**No. 438HT**
- Hand Grinder with Worm Feed..... **No. 438HW**
- Tooth Feed only**No. 438TF**
- Worm Feed only..... **No. 438WF**