

Seedburo Digital Fractional Caliper

These rugged stainless steel calipers have a tri-mode digital display that measures in decimal inches, millimeters, and fractional inches and are ideal for accurate seed measurement or quick screen perforation verification. The units feature a free-wheeling friction roller that advances and retracts the slide bar, scale is graduated in inches and millimeters, auto on/off and includes an Alkaline LR44/257 battery. Accuracy: ±0.001" (±0.02mm). Resolution and Repeatability: 0.0005" or 0.01mm or 1/64". Available in 3 sizes: 6" (153mm), 8" (203mm), 12" (305mm).



147-6

- Digital Fractional Caliper, 6"No. 147-6
- Digital Fractional Caliper, 8"No. 147-8
- Digital Fractional Caliper, 12"No. 147-12

Seedburo Lighted Tweezers

Shed some light on your seeds with these Lighted Tweezers. The tweezers feature a push-button on/off powerful LED light, non-slip comfort grips, heat-treated steel tooling and corrosive resistant non-glare black finish.



Left to right: 70401, 70403, 70408

- Lighted Tweezers, 6 1/4" (159mm) Smooth Sharp PointNo. 70401
- Lighted Tweezers, 6 1/4" (159mm) Serrated Blunt PointNo. 70403
- Lighted Tweezers, 6 1/4" (159mm) Serrated Medium Point with Curved TipNo. 70408

Seedburo Pocket Lighted Magnifier

This powerful magnifier with optional light (via on/off switch) is ideal for today's seed analyst or grain inspector. The unit features an unbreakable, scratch resistant 1 1/2" lens with 7x magnification, an LED light source and fits comfortably in your pocket. The magnifier light requires 3-AAA batteries (not included). Ship wt: 1 lb, Ship dims: 7" (L) x 5" (W) x 3" (H).



- Pocket Lighted 7x MagnifierNo. PLM7X

Seedburo Tweezer Magnifier

This tweezer magnifier is useful for detailed projects. The unit has a 5-power lens, stainless steel tweezer and double hinged adjustable arm. Made with a durable plastic handle and frame with a fine tip. Tweezer length 2 1/4".



- Tweezer MagnifierNo. TWEZEMAGNIFIER

Seedburo Lighted Magnifier

This lighted magnifier puts the light directly on your subject to be viewed. The unit has two magnifications, 3x and 6x, as well as an integrated stand for hands-free use. Push on/Push off powerful LED light, contoured handle for handheld use. 3 1/2" diameter lens. 3x magnification (6x inset lens). Operated by 2-AA batteries included. Includes protective pouch.



- Lighted MagnifierNo. LIGHTMAGNIFIER