

Special Perforation Sieves

Seedburo offers these special perforation sieves with cold-rolled 20ga steel perforated discs spun into standard aluminum sides. We offer the following perforations: round, oblong (rounded ends), slotted (square ends), and herringbone patterns. In addition, we have some perforations available in millimeters. Please request other sizes that may not be listed. These dockage sieves are recommended for use with the Seedburo Sieve Shaker found on [page 40](#). **MADE IN USA.**



Round

| Perforation Size | Part No. | Perforation Size | Part No. | Perforation Size | Part No. | Perforation Size | Part No. |
|---------------------|----------|----------------------|----------|---------------------|----------|----------------------|----------|
| 1/25" Round..... | No. 0022 | 9.5/64" Round | No. 0038 | 17.5/64" Round..... | No. 0062 | 28.5/64" Round | No. 0071 |
| 1/24" Round | No. 0021 | 10/64" Round..... | No. 0036 | 18/64" Round..... | No. 0065 | 29/64" Round | No. 0073 |
| 1/23" Round..... | No. 0024 | 10.5/64" Round..... | No. 0039 | 18.5/64" Round..... | No. 0066 | 30/64" Round | No. 0093 |
| 1/22" Round..... | No. 0019 | 11/64" Round | No. 0040 | 19/64" Round..... | No. 0067 | 32/64" Round | No. 0092 |
| 1/21" Round | No. 0023 | 11.5/64" Round | No. 0041 | 19.5/64" Round..... | No. 0068 | 33/64" Round | No. 0078 |
| 1/20" Round..... | No. 0020 | 12/64" Round..... | No. 0042 | 20/64" Round | No. 0070 | 34/64" Round | No. 0086 |
| 1/18" Round..... | No. 0018 | 12.5/64" Round..... | No. 0047 | 21/64" Round..... | No. 0072 | 35/64" Round | No. 0079 |
| 1/15" Round..... | No. 0015 | 13/64" Round..... | No. 0043 | 22/64" Round | No. 0075 | 36/64" Round | No. 0081 |
| 1/14" Round | No. 0010 | 13.5/64" Round..... | No. 0048 | 22.5/64" Round..... | No. 0076 | 38/64" Round | No. 0087 |
| 1/13" Round..... | No. 0005 | 14/64" Round..... | No. 0045 | 23/64" Round | No. 0077 | 40/64" Round | No. 0094 |
| 4/64" Round | No. 0025 | 14.5/64" Round..... | No. 0046 | 24/64" Round | No. 0080 | 42/64" Round | No. 0084 |
| 7/64" Round..... | No. 0035 | 15/64" Round..... | No. 0050 | 25/64" Round | No. 0085 | 48/64" Round | No. 0083 |
| 7.5/64" Round..... | No. 0037 | 15.5/64" Round..... | No. 0051 | 25.5/64" Round..... | No. 0069 | 56/64" Round | No. 0096 |
| 7.75/64" Round..... | No. 0033 | 16/64" Round..... | No. 0055 | 26/64" Round | No. 0090 | 1" Round | No. 0101 |
| 8.5/64" Round | No. 0031 | 16.5/65" Round..... | No. 0063 | 27/64" Round..... | No. 0091 | | |
| 9/64" Round | No. 0032 | 17/64" Round..... | No. 0060 | 28/64" Round | No. 0095 | | |

Herringbone, Oblong and Slotted

| Perforation Size | Part No. | Perforation Size | Part No. | Perforation Size | Part No. | Perforation Size | Part No. |
|-----------------------------|----------|-----------------------------|----------|------------------------------|------------|----------------------------|----------|
| 1/24" x 1/2" Oblong..... | No. 0124 | .076" x 3/4" Oblong | No. 0145 | 10/64" x 3/4" Oblong..... | No. 0141 | 14/64" x 3/4" Oblong..... | No. 0315 |
| 1/22" x 1/2" Oblong..... | No. 0122 | 5.5/64" x 3/4" Oblong..... | No. 0104 | 10.5/64" x 3/4" Oblong | No. 0108 | 15/64" x 3/4" Oblong..... | No. 0117 |
| 1/20" x 1/2" Oblong..... | No. 0128 | 5.75/64" x 3/4" Oblong..... | No. 0243 | 11/64" x 3/4" Oblong..... | No. 0110 | 16/64" x 3/4" Oblong..... | No. 0118 |
| 3/64" x 5/16" Oblong..... | No. 0121 | 6/64" x 3/4" Oblong | No. 0088 | 11/64" x 3/4" HB..... | No. 0110HB | 17/64" x 3/4" Oblong | No. 0129 |
| 3.5/64" x 1/2" Oblong..... | No. 0127 | 6.5/64" x 3/4" Oblong | No. 0099 | 11.5/64" x 3/4" Oblong..... | No. 0111 | 18/64" x 3/4" Oblong..... | No. 0112 |
| 1/18" x 3/4" Oblong..... | No. 0142 | 7/64" x 3/4" Oblong | No. 0097 | 12/64" x 3/4" Oblong..... | No. 0115 | 19/64" x 3/4" Oblong..... | No. 0119 |
| 1/16" x 3/4" Oblong..... | No. 0143 | 7.5/64" x 3/4" Oblong..... | No. 0098 | 12/64" x 3/4" HB..... | No. 0115HB | 20/64" x 3/4" Oblong..... | No. 0125 |
| 1/14" x 1/2" Oblong..... | No. 0132 | 8/64" x 3/4" Oblong | No. 0100 | 12.5/64" x 3/4" Oblong | No. 0120 | 21/64" x 3/4" Oblong..... | No. 0146 |
| 4/64" x 1/2" Oblong | No. 0109 | 8.5/64" x 3/4" Oblong | No. 0114 | 13/64" x 3/4" Oblong..... | No. 0116 | 22/64" x 3/4" Oblong..... | No. 0138 |
| 4.5/64" x 3/4" Oblong..... | No. 0137 | 9/64" x 3/4" Oblong | No. 0105 | 13/64" x 3/4" HB..... | No. 0116HB | 26/64" x 3/4" Oblong..... | No. 0147 |
| 4-7/8/64" x 3/4" Oblong.... | No. 0113 | 9.5/64" x 3/4" Oblong | No. 0107 | 13.5/64" x 3/4" Oblong | No. 0126 | 32/64" x 3/4" Oblong..... | No. 0139 |

Millimeter

| Perforation Size | Part No. |
|------------------|----------|
| 6mm Round..... | No. 0150 |
| 7mm Round..... | No. 0151 |
| 8mm Round..... | No. 0152 |
| 9mm Round..... | No. 0153 |
| 10mm Round | No. 0154 |

Wire Mesh

| Perforation Size | Part No. |
|-------------------------|----------|
| 36 x 36 Wire Mesh | No. 0135 |
| 40 x 40 Wire Mesh | No. 0140 |
| 4 x 18 Wire Mesh..... | No. 0148 |
| 4 x 20 Wire Mesh..... | No. 0144 |

Triangle

| Perforation Size | Part No. |
|---|----------|
| 10/64" Triangle (9/64" Inscribed Circle) | No. 0102 |
| 11/64" Double Triangle (6.5/64" Inscribed Circle) | No. 0134 |
| 12/64" Double Triangle (7/64" Inscribed Circle)..... | No. 0103 |