

Equipoise Compact Balance

Ideal for classrooms, laboratories and industrial applications. These seven-button balances allow quick access to 13 different units of measure, percentage, counting and accumulating. All units include automatic zero tracking, automatic tare and an accumulation facility which allows the count to be stored and recalled as an accumulated total.

A detachable enclosed weighing chamber is included in the capacities below 1000 grams. Equipoise can be used with the included AC adapter or an on-board rechargeable battery that offers up to 40 hours of continuous operation. RS-232 output is standard on each model. RS-232 cable is optional.



165341

Standard Features:

- Capacities between 300–3000g
- LCD display
- Quick response
- 13 different units of measure
- Modes include accumulation, counting and percentage
- Detachable enclosed weighing chamber (in capacities <1kg)
- Battery operation for up to 40 hours of continuous use (without backlight)
- Standard RS-232 interface for connecting to a printer or PC

Model	No. 165338	No. 165339	No. 165340	No. 165341
Capacity (g)	300g	600g	1500g	3000g
Resolution (g)	0.01g	0.02g	0.05g	0.1g
Repeatability (g) (±)	0.01g	0.02g	0.05g	0.1g
Linearity (g) (±)	0.02g	0.04g	0.1g	0.2g
Stabilization Time (seconds)	2			
Weighing Units	Gram, Carat, Ounce Troy, Pound, Pennyweight, Grain, Momme, Hong Kong Tael, Singapore Tael, Taiwan, Tael, Tola			
Operating Modes	Weighing, Percent Weighing, Parts Counting, Accumulation			
Pan Size	4.7" (120mm)		4.8" x 5.6" (122mm x 142mm)	
Ship Weight and Dimensions	7 lb (3kg), 7.1" x 8.7" x 3.3" (180mm x 220mm x 85mm)			
Power	AC adapter (12v/500mA), internal rechargeable battery (6V/1.2Ah)			
Warranty	One-year limited warranty			

Ohaus Scout STX™ Touchscreen Scales

Ideal for laboratory and industrial applications, the Ohaus Scout features an easy-to-navigate color touchscreen display and a slim design which allows stackable storage. Standard feature includes: quick and efficient navigation with informative color touchscreen display. Fast weighing speed and high resolution with stabilization time as fast as 1 second. Integrated superior overload protection system rated at 10 times the capacity. Multiple connectivity options including RS232, USB Host, USB Device, Ethernet and Bluetooth (available as an accessory). AC adapter included or powered by 4-AA batteries (not included). Ship wt: 4 lbs, Ship dims: 12" (L) x 10" (W) x 6" (H).



Scout STX

Model	No. STX223*	No. STX422	No. STX621	No. STX8200
Capacity	220g	420g	620g	8200g
Readability	0.001	0.01	0.1	1
Pan Size	0.37" (93mm)	0.47" (120mm)	6.7" x 5.5" (170 x 140mm)	6.7" x 5.5" (170 x 140mm)

*Requires draftshield