



Triple Beam Balance

Seedburo offers the very best in mechanical balances. Manufactured by Ohaus, these balances carry a 5 year warranty. All of these scales are ideal for many non-legal for trade applications. The triple beam balances feature a magnetic damper for fast, accurate weighing and easy manipulation of the poises.

The capacity of the gram scale is 610g. With the addition of the **No. 707** weight set, it increases the capacity to 2610g. The front beam reads 10g x 0.1g, rear beam 100g x 10g and the center beam is 500g x 100g. The mechanical series is available with a poly (**No. 99**) or stainless steel scoop (**No. 99S**) which measures 12" x 6" x 2 3/4", or with round 6" diameter stainless steel plate (**No. 760**). Ship wt: 6 lbs, Ship dims: 21" (L) x 9" (W) x 10" (H).

- Gram Scale with Poly Scoop, Spider Platform Style **No. 99**
- Gram Scale with Stainless Steel Scoop, Spider Platform Style **No. 99S**
- Gram Scale with Stainless Steel Plate and Tare Beam **No. 760**



Dial-O-Gram Balances

The Dial-O-Gram Balances eliminate the front beam of a Triple Beam Balance and replace it with a convenient direct reading dial that reads 10g x 0.1g. These scales are equipped with magnetic damping and are suited for non-legal for trade applications.

Capacity is 610g and 2610g with a **No. 707** weight set. The stainless plate is 6" in diameter to accommodate the footed stainless scoop and counter weight (**No. 703S**) Dial 10g x 0.1g, rear beam 100g x 10g, center beam 500g x 100g. Ship wt: 9 lbs, Ship dims: 22" (L) x 9" (W) x 9" (H).

- Dial-O-Gram with Stainless Steel Plate and Tare Beam **No. 1650-00**
- Dial-O-Gram with Stainless Steel Plate, Tare Beam and No. 707 Weight Set **No. 1650-W0**

Triple Beam and Dial-O-Gram Accessories



- Weight Set, Two 1000G and One 500G Weights..... **No. 707**
- Stainless Flat Bottom Scoop and Counter Weight..... **No. 703S**
- Poly Scoop and Counter Weight **No. 703**
- Footed Poly Scoop Used with No. 99 Gram Scale (White) **No. 1101**
- Footed Poly Scoop Used with No. 99 Gram Scale (Gray) **No. 1101G**
- Stainless Steel Scoop Used with No. 99S Gram Scale..... **No. 80250400**
- Maxi-Scoop, 10" x 2 1/2" Polypropylene, for Tare Models Only **No. 3369**



Ohaus Defender™ D500M Portable Platform Scale

THE OHAUS Defender™ D500M is a portable, mechanical industrial scale designed to meet a variety of applications where simplicity and economy are a must. The D500M features an easy-to-read, dual-side flat beam with a shadow-proof poise, providing quick reading under most lighting conditions. The scale is constructed with die-cast beam, fabricated carbon steel base with 4" x 5" diameter phenolic wheels. The platform is 28.13" x 19.25" and the scale capacity is 1000 lbs x 0.5 lb. The D500M is NTEP Certified. Net wt: 114 lbs, Ship wt: 124 lbs, Ship dims: 50" (L) x 32" (W) x 8" (H).

- D500M Portable Platform Scale..... **No. D500M**