



**Seedburo Digital Psychrometer with Dew Point and Wet Bulb**

This certifiable digital psychrometer offers a solution to a wide range of applications in monitoring environmental conditions. Displays wet bulb, dew point, humidity and temperature on a large digital display. No twirling or charts necessary. Includes Min/Max memory with reset key. Range: -4° to 122° F (-30° to 50°C). Temperature accuracy: ±1.8°F. Humidity accuracy: ± 5% reading. Unit runs on two AAA alkaline batteries (included). Net wt: 1 lb, 6.7" (L) x 1.85" (W) x 0.5" (H), Ship wt: 2 lbs, Ship dims: 12" (L) x 10" (W) x 3" (H).

Digital Psychrometer with Dew Point and Wet Bulb .....No. SAM990DW



**Seedburo Calibratable 7 Function Digital Psychrometer**

This combination meter measures and displays ambient temperature, relative humidity (RH), both dew point and wet bulb temperatures on a bright dual-readout LCD display. Includes an external stem probe for measuring surface or internal temperature of a solid object or the temperature of a liquid or gas. NIST calibration available. F°/C° switch and an adjustable auto power off trigger. Unit includes one 9V alkaline battery. Net wt: 3.7 oz, Ship wt: 2 lbs, Actual dims: 7.56" (L) x 2.4" (W) x 1.57" (H), Ship dims: 12" (L) x 10" (W) x 3" (H).

Calibratable 7 Function Psychrometer .....No. EP8710

Calibration Reference 33% and 75% Salts for Field Calibration .....No. HR33 or HR75



**Seedburo Sling Psychrometer**

The Seedburo Sling Psychrometer determines relative humidity by measuring evaporation of water into the surrounding air. The unit compares wet and dry bulb temperatures. The %RH is read from the slide rule calculator integrated into the psychrometer. Standard humidity tables can also be used. Unit contains two spirit filled thermometers and wick. Net wt: 4 oz, Actual dims: 8½" (L) x 1½" (W) x ⅝" (H), Ship wt: 1 lb, Ship dims: 12" (L) x 3" (W) x 10" (H).

Sling Psychrometer, 20 to 120°F.....No. 318

Sling Psychrometer, -5 to 50°C.....No. 318C

Replacement Thermometer °F.....No. 318T

Replacement Thermometer °C .....No. 318TC

Replacement Wicks.....No. 318W



**Delmhorst Thermo-Hygrometer**

The Delmhorst Thermo-Hygrometer, offered by Seedburo Equipment Company, is an effective way to control the effects of humidity and temperature in a variety of locations such as seed storage rooms, walk in germinators or coolers, or grain drying rooms. The HT-3000 displays relative humidity from 0–100% with a ± 3% accuracy and temperature from -4° F to 158° F (-20° C to 70° C) with up to ±1.8° accuracy. Dew point range is -40° F to 158° F (-40° C to 70° C). The HT-3000 has a quick response time, selectable display screens, data hold feature, graphic display with backlight, auto power down when the unit has been idle and hard rugged case. Actual wt: 1 lb., Actual dims: 3" (L) x 1½" (W) x 8" (H), Ship wt: 3 lbs, Ship dims: 13" (L) x 10" (W) x 7" (H).

Thermo-Hygrometer .....No. HT-3000



**Seedburo Mason Hygrometer**

The Mason Hygrometer has easy reading thermometers for wet and dry bulb temperatures. These spirit-filled thermometers have a dual temperature range of 20 to 120°F and -5 to 50°C and are mounted on a durable plastic housing for convenient wall mounting. Unit includes humidity tables for easy reference to atmospheric humidity. Net wt: 8 oz, Actual dims: 3" (W) x 12" (H), Ship wt: 1 lb, Ship dims: 12" (L) x 3" (W) x 10" (H).

Mason Hygrometer .....No. 312D