

**Achieva Precision Table Top Germinator**



Seedbuero's Achieva tabletop germinator offers features similar to our larger units, but these units will save you both space and money! The ATTGPT is ideal for small volume seed labs that require precision, but do not require a full size germinator.

The ATTGPT germinator is equipped with a programmable computer touch screen controller. Large full color display for easy programming and display of: test and cycle(s), time remaining, program recall, editing, time of day, test duration, run/program or auto/manual modes, real time clock, data download function, heat/cool and light/dark levels. Lights can be on or off in either mode. A digital display data logger stores both temperature and humidity readings at the desired time interval and duration. Data is easily downloaded to any Windows®-based PC or laptop computer. Software and download cable included. Each unit is equipped with interior cool white bright LED light strips to provide for complete light/dark (lights are not adjustable) at 180 lumens, 6500k. Interior has reflective surfaces for greater light spread.

Temperature range is between 5°C–40°C, ± 1.1°C (41°F–104°F, ± 2°F). A Digital temperature display shows both set-point and actual cabinet temperature. An internal fan distributes air evenly and gently with minimal temperature stratification. Heat is maintained by a long-life, stainless steel, and low watt density heating element. Uniform cooling is maintained through the use of an air cooled hermetically sealed condenser compressor unit. The ATTGPT is a dry style germinator designed for containerized germination testing. Two open wire shelves, 20" x 14" at 2½" spacing, are included. Interior dims: 19" (L) x 20.75" (W) x 24.5" (H). Net wt: 75 lbs, Actual dims: 20" (L) x 23.14" (W) x 39.53" (H), Ship/Export wt: 285 lbs, Ship/Export dims: 36" (L) x 41" (W) x 50" (H).

|   |               |
|---|---------------|
| Achieva Table Top Germinator, 115V, 60Hz.....             | No. ATTGPT/B  |
| Achieva Table Top Germinator, 220V, 50Hz.....             | No. ATTGPT/C  |
| Crating for Domestic Shipping or Exporting, Required..... | No. ATTGPT/EX |



**Seedbuero Light-Dark Germinator**

A clear, angled plexiglass top allows light in and any condensation formed drains off without dripping on the samples. Stainless steel inside and out. High-density foam provides a non-absorbent, insect free thermal barrier to provide an even temperature within the chamber. Eight aluminum 15" x 21" trays are included, providing a total of 2520 square inches of testing area. A water/ice pan is included. Fill with water or ice when additional humidity or below ambient temperatures are required. Heat is set by a digital controller and can be set from 5°C–40°C (41°F–104°F). Digital display and rubber feet are standard. Interior dims: 16.9" (D) x 18.65" (H) x 21.82" (W), Net wt: 90 lbs, Actual dims: 24.08" (L) x 19.87" (W) x 26.51" (H), Ship/Export wt: 320 lbs, Ship/Export dims: 36" (L) x 32" (W) x 54" (H).

|   |              |
|---|--------------|
| L-D Germinator, 110V, 50/60Hz, 1PH..... | No. 548      |
| L-D Germinator, 220V, 50/60Hz, 1PH..... | No. 548/A    |
| Replacement Heating Element, 110V.....  | No. 54H      |
| Replacement Heating Element, 220V.....  | No. 54H/A    |
| Replacement Thermostat.....             | No. 54TR     |
| Replacement Digital Controller.....     | No. 54DC     |
| Replacement Thermometer.....            | No. 548T     |
| Replacement Tray.....                   | No. 548-TRAY |