



★ **MADE IN USA**
757 with 757-M tube set, sold separately

Seedbuero South Dakota Seed Blower

The 757 South Dakota Seed Blower utilizes airflow through a column generated by a blower motor combination. Separations are accurately controlled by a calibrated valve cap in the top of the column. We offer three column sizes: small 1½" I.D. x 23¼" H, medium 3" I.D. x 34¼" H, and large 4" I.D. x 35¼" H. Air columns are precision built to trap lightweight seed and chaff, collecting pure seed in the "cup" located at the bottom of the column. The cup has a wire mesh bottom to collect the seed and allow the air to flow through the seed. Mesh sizes are: small (100 mesh), medium (50 mesh) and large (30 mesh). The 757 base measures 16" (L) x 21" (W) x 32" (H). This permits easy viewing of the calibrated air valves on all three-column sizes. Additional features include a 5-minute analog timer, on/off control knob and a hand brake to stop the motor instantly to eliminate coasting of the blower fan. Net wt: 110 lbs, Ship wt: 140 lbs, Ship dims: 25" (L) x 21" (W) x 41" (H). **Tube Sets:** Ship wt: 6 lbs, Ship dims: 21" (L) x 9" (W) x 9" (H).

Seedblower Table, 120V, 50/60Hz, 3450 RPM, 1/3HP Motor.....	No. 757
Seedblower Table, 220V, 50/60Hz, 3450 RPM, 1/3HP Motor.....	No. 757/A
Small Tube Set (1½" column).....	No. 757-S
Medium Tube Set (3" column).....	No. 757-M
Large Tube Set (4" column)	No. 757-L
Required Crating	No. 757/DOM



★ **MADE IN USA**
757-RC

Seedbuero South Dakota Seed Blower Cup Set

Also available is the 757-RC, 4" tube with removable cup for easy access to the cleaned sample. It has removable containers off the side of the column, in lieu of baffles as in the standard columns. Net wt: 8 lbs, Ship wt: 10 lbs, Ship dims: 21" (L) x 21" (W) x 16" (H).

4" Tube with Removable Cup 30 Mesh	No. 757-RC
Extra Cup with 100 Mesh for 757-RC Tube Set	No. 757-CUP100
Extra Cup with 50 Mesh for 757-RC Tube Set.....	No. 757-CUP50
Extra Cup with 30 Mesh for 757-RC Tube Set.....	No. 757-CUP30



★ **MADE IN USA**
SABSC/B

Air Blast Seed Cleaner

This cleaner is a tabletop model for final cleaning of small samples. Grain first passes through a screening tray for removal of large debris, then into an air chamber where final separation takes place. Clean grain is collected in a removable catch pan. Chaff and light debris are blown into a second catch pan for inspection and re-cleaning. Standard equipment includes air-intake gate on blower and 2 seed catch pans. Screen sizes ½", ¼", ⅛", ⅜" (not included) wood base, 6 ft. drop cord and switch. Ship wt: 98 lbs, Ship dims: 24" (L) x 18" (W) x 32" (H). Export wt: 100 lbs, Export dims: 25" (L) x 24" (W) x 32" (H). **Screen:** Ship wt: 7 lbs, Ship dims: 18" (L) x 12" (W) x 12" (H).

Seed Cleaner 115V, 60Hz without Screens.....	No. SABSC/B	Screen with ⅜" Opening.....	No. SABSC-S38
Seed Cleaner 220V, 50Hz without Screens.....	No. SABSC/C	Seed Catch Pans.....	No. SABSCP
Screen with ½" Opening.....	No. SABSC-S12	Domestic Crating.....	No. SABSC/DOM
Screen with ¼" Opening.....	No. SABSC-S14	Export Crating	No. SABSC/EX
Screen with ⅛" Opening	No. SABSC-S18		