

Hand Testing Screens

Seedbuero's hand testing screens are a must for anyone who buys, cleans, or grades seed or grain. They are made of selected straight grained wood with nesting frames. These screens can be stacked together with a blank bottom pan to make several separations at one time.

A generous assortment of sizes will enable you to experiment with various separations and select the correct screens for your larger cleaning or sizing

operations. Our hand test screens are an excellent choice for experimental or laboratory use when only small samples are available. Net wt: 1/2 lb, Ship wt: 3 lbs. All hand testing screens are **MADE IN USA**.

- Hand Test Screens, 12" x 12".....**No. 98-12Q**
- Hand Test Bottom Pan, 12" x 12".....**No. 98-12BQ**
- Screen Perforation Gauge.....**No. PFG**

*** Only available in metric size*

Round Holes (64ths of an inch)		Round Holes (64ths of an inch)		Oblong Holes		Oblong Holes		Wires/Sq. Inch
Inches	MM	Inches	MM	Inches	MM	Inches	MM	10 x 10
1/25	1.016	23	9.127	1/26 x 1/2	0.98 x 12.70	13/64 x 3/4	5.159 x 19.050	12 x 12
1/24	1.058	23 1/2	9.330	1/25 x 1/2	1.02 x 12.70	13/64 x 1 1/4	5.159 x 31.75	14 x 14
1/23	1.102	24	9.525	1/24 x 1/2	1.06 x 12.70	13 1/2/64 x 3/4	5.357 x 19.050	15 x 15
1/22	1.153	24 1/2	9.720	1/23 x 1/2	1.10 x 12.70	13 1/2/64 x 1 1/4	5.357 x 31.75	16 x 16
3/64	1.190	25	9.922	1/22 x 1/2	1.153 x 12.70	14/64 x 3/4	5.556 x 19.050	18 x 18
1/21	1.200	25 1/2	10.120	3/64 x 5/16	1.191 x 7.938	14 1/2/64 x 3/4	5.75 x 19.050	20 x 20
1/20	1.270	26	10.319	3/64 x 1/2	1.19 x 12.70	15/64 x 3/4	5.953 x 19.050	22 x 22
1/19	1.336	26 1/2	10.520	1/20 x 1/2	1.270 x 12.700	15 1/2/64 x 3/4	6.15 x 19.050	24 x 24
1/18	1.410	27	10.716	1/19 x 1/2	1.34 x 12.70	16/64 x 3/4	6.350 x 19.050	26 x 26
1/16	1.587	27 1/2	10.910	1/18 x 3/4	1.410 x 19.050	16 1/2/64 x 3/4	6.55 x 19.050	28 x 28
1/15	1.691	28	11.113	1/18 x 1/2	1.410 x 12.700	17/64 x 3/4	6.747 x 19.050	30 x 30
1/14	1.813	28 1/2	11.310	1/17 x 1/2	1.494 x 12.700	17 1/2/64 x 3/4	6.95 x 19.050	32 x 32
1/13	1.953	29	11.509	1/15 x 1/2	1.691 x 12.700	18/64 x 3/4	7.144 x 19.050	34 x 34
1/12	2.116	29 1/2	11.710	1/15 x 3/4	1.69 x 19.05	18 1/2/64 x 3/4	7.34 x 19.050	36 x 36
6	2.381	30	11.906	4 1/2/64 x 1/2	1.785 x 12.700	19/64 x 3/4	7.541 x 19.050	38 x 38
6 1/2	2.579	30 1/2	12.100	4 1/2/64 x 3/4	1.785 x 19.050	20/64 x 3/4	7.939 x 19.050	40 x 40
7	2.778	31	12.303	4 3/4 x 3/4	1.890 x 19.050	21/64 x 3/4	8.33 x 19.050	45 x 45
7 1/2	2.976	31 1/2	12.500	4 7/8/64 x 1/2	1.934 x 19.050	22/64 x 3/4	8.731 x 19.050	50 x 50
7 3/4	3.080	32	12.700	4 7/8/64 x 3/4	1.934 x 19.050	22 1/2/64 x 1	8.93 x 25.400	55 x 55
8	3.175	33	13.100	5/64 x 1/2	1.984 x 12.700	23/64 x 3/4	9.13 x 19.050	60 x 60
8 1/2	3.373	34	13.949	5/64 x 3/4	1.984 x 19.050	24/64 x 3/4	9.525 x 19.050	3 x 20
8 3/4	3.470	34 1/2	13.690	5.25/64 x 3/4	2.080 x 19.050	28/64 x 3/4	11.110 x 19.050	4 x 18
9	3.572	35	13.890	5 1/2/64 x 1/2	2.182 x 12.700	30/64 x 3/4	11.910 x 19.050	4 x 19
9 1/2	3.770	35 1/2	14.090	5 1/2/64 x 3/4	2.182 x 19.050	32/64 x 3/4	12.700 x 19.050	4 x 20
9 3/4	3.870	36	14.287	5.8/64 x 3/4	2.300 x 19.050	36/64 x 3/4	14.29 x 19.050	4 x 22
10	3.969	38	15.081	6/64 x 1/2	2.380 x 12.700	**	2.15mm x 20mm	6 x 18
10 1/2	4.167	40	15.875	6/64 x 3/4	2.380 x 19.050	**	2.25mm x 20mm	6 x 20
10 3/4	4.270	42	16.669	6/64 x 1 1/32	2.380 x 26.160	**	3.75mm x 20mm	6 x 22
11	4.366	44	17.460	6 1/2/64 x 3/4	2.579 x 19.050	**	4mm x 20mm	6 x 24
11 1/2	4.564	46	18.260	7/64 x 3/4	2.778 x 19.050	**	4mm x 23mm	6 x 26
11 3/4	4.660	47	18.650	7.25 x 3/4	2.88 x 19.050	**	4mm x 31mm	6 x 28
12	4.763	48	19.050	7 1/2/64 x 3/4	2.972 x 19.050	**	4.25mm x 20mm	6 x 30
12 1/2	4.961	50	19.840	8/64 x 3/4	3.175 x 19.050	**	4.25mm x 31mm	6 x 32
13	5.159	52	20.640	8/64 x 1	3.18 x 25.4	**	4.35mm x 20mm	6 x 34
13 1/2	5.357	54	21.430	8/64 x 1 1/4	3.18 x 31.75	**	4.35mm x 31mm	6 x 36
14	5.556	55	21.83	8 1/2/64 x 3/4	3.373 x 19.050	**	4.4mm x 20mm	6 x 38
14 1/2	5.754	56	22.225	8 1/2/64 x 1 1/4	3.37 x 31.75	**	4.4mm x 31mm	6 x 40
14 3/4	5.850	57	22.620	9/64 x 1/2	3.57 x 31.75	**	4.5mm x 20mm	6 x 42
15	5.953	64	25.400	9/64 x 3/4	3.527 x 19.050	**	4.5mm x 31mm	
15 1/2	6.150	1 3/4	31.750	9/64 x 1 1/4	3.57 x 31.75	**	4.75mm x 20mm	
16	6.350	**	5.5mm	9 1/2/64 x 3/4	3.770 x 19.050	**	4.75mm x 31mm	Triangular Holes
16 1/2	6.548	**	7.5mm	9 1/2/64 x 1 1/4	3.77 x 31.75	**	5.25mm x 20mm	9/64 or 5 1/2V
17	6.747	**	10mm	10/64 x 3/4	3.969 x 19.050	**	5.25mm x 31mm	10/64 or 6V
17 1/2	6.945			10/64 x 1 1/4	3.969 x 31.75	**	5.5mm x 20mm	11/64 or 6 1/2V
18	7.144			10 1/2/64 x 3/4	4.17 x 19.050	**	5.5mm x 31mm	12/64 or 7V
18 1/2	7.342			10 1/2/64 x 1 1/4	4.17 x 31.75	**	5.75mm x 20mm	17/64 or 9 1/2V
19	7.541	32/64	12.700	11/64 x 3/4	4.366 x 19.050	**	5.75mm x 31mm	
19 1/2	7.739			11/64 x 1 1/4	4.366 x 31.75	**	6mm x 20mm	
20	7.938			11 1/2/64 x 3/4	4.564 x 19.050	**	6mm x 23mm	
20 1/2	8.136			11 1/2/64 x 1 1/4	4.564 x 31.75	**	6mm x 31mm	
21	8.334			12/64 x 3/4	4.763 x 19.050	**	12mm x 30mm	
21 1/2	8.532			12/64 x 1 1/4	4.763 x 31.75			
22	8.731			12 1/2 x 3/4	4.961 x 19.050			
22 1/2	8.929			12 1/2 x 1 1/4	4.961 x 31.75			