



★ MADE IN
★ USA

400/B with optional
Back Air Catchall with removable trays

Clipper Office Tester and Cleaner

The Clipper Model 400 Office Tester and Cleaner offered by Seedbuero is an extremely efficient cleaner for office and laboratory use. This machine will provide accurate, dependable performance for small lot cleaning, testing or sampling. Approximate capacity would be 1 bushel per hour for most crop seeds. It is commonly used in seed plants for sampling large lots of seed for correct screen selection on larger, commercial size cleaners. It is also widely used in university laboratories for agronomy research and testing. Flower and vegetable growers use the Model 400 for specialized cleaning, grading and sizing. The Model 400 has a compact, portable design and features an eccentric (shoe) drive that uses two (2) 10" x 15" screens. One top screen for scalping, one bottom screen for cleaning and grading. Perforated metal or wire cloth screens

are available in over 175 different sizes (For screen perforations, see **page 86**). Other features include a powerful bottom blast fan driven by a three-step pulley, adjustable slide on the air intake openings to provide intermediate air regulation and an optional, removable catchall with either removable trays (shown) or solid bottom to catch air screenings. Net wt: 70 lbs, Actual dims: 25" (L) (45¾" with catchall) x 17" (W) x 23" (H), Ship wt: 110 lbs, Ship dims: 32" (L) x 24" (W) x 32" (H), Export wt: 185 lbs, Export dims: 33" (L) x 25" (W) x 33" (H).

- Office Tester and Cleaner Complete with 2 Screens **No. 400/B**
1/6HP, 115V, 60Hz/1Ph
- Office Tester and Cleaner Complete with 2 Screens **No. 400/E**
1/6HP, 230V, 60Hz/1Ph
- Office Tester and Cleaner Complete with 2 Screens **No. 400/C**
1/4HP, 230V, 50Hz/1Ph
- Office Tester and Cleaner Complete with 2 Screens **No. 400/CE**
1/4HP, 230V, 50Hz/1Ph Motor, CE Certified

- Back Air Catchall with Solid Bottom (Optional) **No. 400BAC**
- Back Air Catchall with Removable Trays (Optional) **No. 400BACT**
- Extra Perforated Metal Screens **No. 400SM**
- Extra Wire Mesh Screens **No. 400SW**
- Required Domestic Packing **No. 400/DOM**
- Required Export Crating **No. 400/EX**



★ MADE IN
★ USA

AGM224

Clipper Model AGM224 Seed and Grain Cleaner

The Clipper AGM224 cleaner is a small cleaner that will provide excellent results and fits the market for those in need of a cleaner between the sample lot size of the Office Clipper and the mid-size capacity of the Eclipse 324 cleaner. This unit is ideal for small production, in lab, or on-farm applications. The AGM224 is a durable cleaner and features a steel frame, galvanized steel shoe, scotch-ply hangers, heavy-duty eccentric (shoe) drive, ball tray cleaning system and plain inlet hopper. Designed with two (2) 24½" x 22¼" screens. The top screen handles scalping, while the bottom screen offers thorough cleaning and grading. An adjustable air separation efficiently completes the cleaning process. With over 175 different screen sizes (perforated metal and wire mesh) available the AGM224 is an extremely versatile machine. Net wt: 575 lbs, Actual dims: 54" (L) x 41" (W) x 67" (H), Ship wt: 895 lbs, Ship dims: 56" (L) x 45" (W) x 70" (H). (For screen perforations, see **page 86**).

Approximate Capacities (Bushels Per Hour)

Commodity	Fast Clean	Precision Clean
Soybeans, Wheat, Corn	*50 BPH	*30 BPH
Rice, Oats, Barley	*40 BPH	*25 BPH

**Note: Capacities may vary depending on seed condition, moisture content, seed varieties, types, and volume of contaminants to be separated and the percentage of foreign material acceptable in the final product.*

- Clipper AGM224 Cleaner Complete with 2 Screens **No. AGM224**
¾HP, 110/220V, 60Hz, 1 Ph
- Clipper AGM224 Cleaner Complete with 2 Screens **No. AGM224/A**
¾HP, 110/220V, 50Hz, 1 Ph
- Extra Standard Metal Perforated Screens **No. AGM224SM**
- Extra Standard Metal Perforated Screens, Cross Slots... **No. AGM224SMCS**

- Extra Square Wire Mesh Screens **No. AGM224SWS**
- Extra Rectangle Wire Mesh Screens **No. AGM224SWR**
- Bagging Auger Attachment **No. AGM224BA**
- Required Domestic Crating **No. AGM224/DOM**
- Required Export Crating **No. AGM224/EX**