

Seedburo Laboratory Oat Huller

The Seedburo Oat Huller is designed to minimize breakage during hulling operation. Hulling a 50 gram sample takes approximately one minute. Hulls and fuzz are separated from the groats and are sorted into two separate collection canisters, allowing for quick weighing of fractions. The fractions allow the lab technician to check for heat and frost damage, mold, weathering, discoloration and odor. The unit's simple and rugged design nearly eliminates the need for spare parts. The unit comes complete with 1/3HP, single phase motor, and requires compressed air at 100 psi. Net wt: 115 lbs, Ship wt: 218 lbs, Ship dims: 27" (L) x 42" (W) x 90" (H). Export wt: 300 lbs, Export dims: 90" (L) x 42" (W) x 28" (H).

Oat Huller, 110V, 60Hz	No. OH5095/B
Oat Huller, 220V, 50Hz	No. OH5095/C
Oat Huller, 110V, 50Hz	No. OH5095/D
Oat Huller, 220V, 60Hz	No. OH5095/E
Required Domestic Crating	No. OH5095/DOM
Required Export Crating	No. OH5095/EX



OH5095/B

Huller/Scarifier

The Huller/Scarifier is designed and constructed to hull/scarify seeds without damage, if operated properly. Upon entering the hulling chamber, the seed is rubbed against a heavy, rubber-lined outer wall by the inside cylinder. This action will hull and/or scarify the toughest of seeds without damage to the seed. Model 2 hulls, scarifies and removes chaff in one operation. Hulling cylinder available in three types: Perf-grip, Linard, Urethane. Huller is available in 50Hz or 60Hz and in 230V, 380V, or 460V. Please advise which voltage and hertz you require. **Model2:** Ship wt: 600 lbs, Ship dims: 50" (L) x 45" (W) x 42" (H), Export wt: 700 lbs, Export dims: 55" (L) x 48" (W) x 45" (H). **Cylinder:** Ship wt: 65 lbs, Ship dims: 26" (L) x 17" (W) x 16" (H).



MODEL2A/F

Features:

- All-Steel Construction
- Settling Chamber and Exhaust Fan
- Variable Speed Drive
- Feed Control Hopper
- Adjustable Air Control

Hulls/Scarifies:

- Clover Seeds
- Grass Seeds
- Quinoa Seeds
- Lespedeza
- Millet
- Alfalfa
- Barley
- Birdfoot Trefoil
- Kudzu
- Any many other seeds

Huller/Scarifier with Aspiration, 5 HP, TEFC 3PH Motor, Please Advise Voltage and Hertz.....	No. MODEL2A
Huller/Scarifier less Aspiration, 5 HP, TEFC 3/60/230, 460 or 575 Motor, Please Advise Voltage and Hertz.....	No. MODEL2
Replacement Cylinder, Perf-grip.....	No. MODEL2P
Replacement Cylinder, Urethane.....	No. MODEL2U
Required Domestic Crating	No. MODEL2/DOM
Required Export Crating	No. MODEL2/EX

Please contact Seedburo for proper cylinder, capacity or free product evaluation.