



Hand-Held Series

Portable Electronic Scales Models HH120D, HH120 and HH320



Convenient, compact and portable weighing capability in the palm of your hand! These little scales with sleek styling and convenient size are ideal to take anywhere—throw it in a brief case or drop it in a backpack. Whether you choose the HH120D model with basic grams weighing or the HH120 and HH320 models with multiple weighing units, the OHAUS Hand-Held series offers outstanding performance and exceptional value.

Standard Features Include:

- Unique hand-held design
- Energy-saving auto shut-off feature
- Protective snap on cover for storage and transport
- Large high contrast display
- Backlit display on HH120 and HH320 models
- Rugged plastic construction

Hand-Held Series Portable Scales

Specifications

Model	HH120D	HH120	HH320
Capacity x Readability (g)	0-60 x 0.1 60-120 x 0.2	120 x 0.1	320 x 0.1
Linearity (g)	0 - 60g: ± 0.1 60 - 120g: ± 0.2	± 0.1	± 0.2
Weighing units	g	g, oz, dwt, oz t	
Zero range	To capacity by subtraction		
Calibration	User calibration		
Stabilization time	3 seconds		
Maximum overload	150% of capacity		
Auto shut-off	5 minutes no activity		
Optimum operating temperature	64°-77°F / 18°-25°C		
Display (in/mm)	LCD (0.4/11)		
Weigh pan size L x W (in/cm)	3.25 x 3.0 / 8.3 x 7.6		
Power requirements	2-AAA alkaline batteries (included)		
Dimensions L x W x H (in/cm)	5.35 x 3.25 x 0.8/13.6 x 8.3 x 2.0		
Net Weight (lb/kg)	0.34 / 0.16		

Standard Features

Low battery indicator, overload and underload indicators, 2-AAA batteries, rugged plastic construction, protective snap-on cover, easy calibration. Button activated luminescent display on HH120 and HH320.

Industry Leading Quality and Support

All OHAUS Hand-Held Scales are manufactured under an ISO 9001:2008 Registered Quality Management System. Our rugged construction and stringent quality control have been hallmarks of all OHAUS products for over a century.

OHAUS CORPORATION

* 7 Campus Drive
Suite 310
Parsippany, NJ 07054 USA

Tel: 800.672.7722
973.377.9000
Fax: 973.944.7177

www.ohaus.com

*With offices throughout Europe,
Asia, and Latin America*

* ISO 9001
Registered Quality Management System



Seedburo Equipment Company • 2293 South Mt. Prospect Road • Des Plaines, IL 60018-1810 USA
Toll-Free 1-800-284-5779 (USA & Canada) • Telephone +1-312-738-3700 / Fax +1-312-738-5329
www.seedburo.com

80774394_A © Copyright OHAUS Corporation

