



THE OHAUS HAND-HELD SERIES

*electronic scales
models HS-120 and HP-320*



Convenient, compact and portable weighing capability in the palm of your hand! These little scales with sleek styling and convenient size are ideal to take anywhere— throw it in a brief case, drop it in a backpack, or slip it in your pocket. Whether you choose the standard HS-120 model with basic grams weighing, or the professional HP-320 with multiple weighing units, the Ohaus Hand-Held Series offers outstanding performance and exceptional value.

- Unique hand-held design
- Energy-saving auto shut-off feature
- Large high contrast display
- Backlit display and protective cover on HP-320
- Rugged plastic construction



Ohaus Hand-Held Series
Electronic Scales

Models HS-120 and HP-320



Model	HS-120	HP-320
Capacity x Readability (g)	0-60 x 0.1 60-120 x 0.2	320 x 0.1
Linearity (g)	0 - 60g: ±0.1 60 - 120g: ±0.2	±0.2
Weighing units	g	g, oz, dwt, oz t
Zero range	To capacity by subtraction	
Calibration	Factory preset	Auto calibration from keypad
Stabilization time	3 seconds	
Maximum overload	150% of capacity	
Auto shut-off	2 minutes no activity	
Optimum operating temperature	64°-77°F/18°-25°C	
Display (in/mm)	LCD (0.4/11)	
Weigh pan size L x W (in/cm)	3.25 x 3.0/8.3 x 7.6	
Power requirements	2-AAA alkaline batteries (included)	
Dimensions L x W x H (in/cm)	5.35 x 3.25 x 0.8/13.6 x 8.3 x 2.0	
Net Weight (lb/kg)	.34/.16	

Standard Features

Low battery indicator, overload and underload indicators, 2-AAA batteries, rugged plastic construction. HP-320 only: Button activated luminescent display, snap on protective molded cover, soft tactile feedback keypad, low battery indicator.

Optional Accessories

Calibration mass for model HP-320, hard shell protective cover for model HS-120

Industry Leading Warranty

All Ohaus hand-held scales are protected by a 12 month warranty on labor and parts. Our rugged construction and stringent quality control have been hallmarks of Ohaus scales and balances since 1907.

Simple To Order

Contact your dealer for product availability.



Seedburo Equipment Company

1022 WEST JACKSON BLVD., CHICAGO, IL 60607-2990
PHONE (312) 738-3700 FAX: (312) 738-5329
TOLL FREE (800) 284-5779
SEE US @ WWW.SEEDBURO.COM

LUS04-99 © Copyright Ohaus Corporation

